Umluftgerät Typ SST-UL V4

Montage- und Bedienungsanleitung /
Installation and Operating Manual /
Instructions de montage et d'utilisation

Ventilatoren

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Technical ventilation

Operating instructions for recirculation fan, type SST-UL V4

Instructions for safe use

The **SST-UL** recirculation fan is used to extract and clean the exhaust air from hazardous material, safety, chemical and gas cylinder cabinets. Ambient temperatures from +10°C to +40°C are permitted; the relative humidity at the place of use must not exceed 55%. The device ventilates the cabinet depending on the required air change rate, filters out hydrocarbons and indicates the current function by means of two lamps each (RED = Error, GREEN = Service). The concentration of hazardous substances in the exhaust air is shown in 25% steps. A concentration of 100% corresponds to the maximum load. The DSA (dynamic saturation adjustment) function corrects the switching point according to the ambient air quality. The device also has two potential-free outputs for these signals. It is suitable for extracting air from cabinets as defined in EN 14470-1 and DIN 12925 Part 2 and conforms to the requirements of TRbF 20, TRGS 510 and GefStoffV. The user is responsible for ensuring safe operation of the system in accordance with TRbF 20 Section 15. The statutory limit values (TRGS 900 "Occupational exposure limit values" and TRGS 903 "Biological limit values") must be observed at the place of use.

PERSONAL INJURY AND MATERIAL DAMAGE

Improper use could result in broken impellers, fatigue failures or fires due to sparks, possibly leading to personal injury or material damage.

Specific hazards

The device conforms to 2014/34/EU (ATEX) (equipment group II, category 3 and temperature classification T4); it is exclusively intended for extracting air from Zone 2 or non-Ex zones and must be installed in a safe area.

Marking: CE  II _/3G c e IIB T4

All forms of use which deviate from or exceed the limits of use described above are considered to be contrary to the intended purpose. The manufacturer / supplier is not liable for any damages resulting from such use; the sole responsibility lies with the user.

No specific hazards for people or the environment arise provided this equipment is used as intended and only with original accessories, and provided this Installation and Operating Manual as well as the information on safe commissioning are observed.

Operating principle, safe installation, commissioning, adjustment, installation / removal, maintenance

The device is shipped ready to use and install with the exception of the control cables, which are not included in the scope of supply because they must be adapted to the specific conditions on site. You must check the electrical specifications (voltage, frequency) and compare them with the information on the rating plate prior to connecting the device to the mains.

The device may only be installed, operated and maintained by appropriately trained, instructed and authorised persons.

- Operating principle: The air containing hazardous substances is discharged from the hazardous material cabinet via a filter system comprised of a coarse filter, a fine filter and an activated carbon filter. 99.9% of all hydrocarbons are retained in this way by adsorption. This micro-controlled device monitors the volume flow and saturation of the activated carbon filter. The operational status is indicated by three LEDs (GREEN, YELLOW and RED) as well as by an audible signal. To meet the requirements of TRbF 20 Annex L, the exhaust function is interrupted as soon as the saturation LED shows "Error". This ensures that no hazardous substances are discharged from the cabinet into the workroom. The interior of the device is automatically flushed in the event of an error in order to reduce the concentration. If flushing does not cause the concentration to fall below the setpoint, the device goes into "alarm". Alarms can only be cancelled by authorised persons. Short concentration peaks (less than 9 seconds) are simply indicated (all saturation lamps on) without the device going into "error". The coarse and fine filters significantly extend the life of the activated carbon filter by retaining dust particles. The actual saturation up to the limit value is indicated by the three yellow LEDs in 25% steps. The quality of the ambient air is shown by another yellow LED. The brightness of this LED varies according to the amount of the hazardous substances.

- Installation and commissioning: The device is attached to the cabinet by means of the connection fitting (accessory). The potential-free signal outputs are connected by the customer if necessary. The device goes into test mode as soon as it is connected to the mains supply with the mains plug. The LEDs, potential-free outputs and sensors are then checked one at a time.
- Function check: A function check should be carried out whenever the device is started up or as needed (see “Troubleshooting”). When the device is connected, or if you press the TEST button, all lamps light up for 3 seconds. Following this, the relays for the potential-free contacts are switched on and off for 3 seconds each. The flow and saturation sensors are checked in the final test sequence. The flow is checked by switching the technical ventilation system on and off (5 seconds on, then 5 seconds off). The LEDs for the flow change colour from RED to GREEN and the result of the test is indicated during the pauses (RED = Error, GREEN = OK). At the end of the test, the device automatically goes into normal mode. The exhaust air monitoring function can also be checked by briefly covering over the exhaust vent in order to throttle the volume flow.
- Adjustment: The exhaust air monitoring function is set to the required flow and saturation at the factory, so that no further adjustments are necessary when you start up the device for the first time.
- Repairs: Repair work may only be carried out by authorised persons.
- Maintenance: The device must be checked by an authorised person as and when necessary (once a year as a minimum) in accordance with TRbF 20 Section 15.4. The filters must be replaced or cleaned at the same time. Replacement filters can be ordered from the supplier. We recommend stocking a set of replacement filters to avoid interrupting operation. The device records the operating time internally and outputs an audible signal (one 3-second beep every hour after 365 days) to indicate that maintenance is due.


Troubleshooting

| Fault | Cause | Remedy |
|--|--|---|
| Device not working | No power Fuse defective Device defective | Check the mains voltage Replace the fuse of the mains connector Send in the device |
| Red flow lamp ON | Dust filter clogged Volume flow interrupted | Replace the coarse filter Replace the fine filter Remove the obstruction from the exhaust vent |
| Red saturation lamp ON, audible signal, air flow interrupted | Activated carbon filter saturated Concentration of hazardous substances too high Other hazardous substances contained in exhaust air | Replace the filter Check the contents of the cabinet and remove any leaked substances Start up the device again (authorised persons only) Contact the supplier |


Specifications

| | | |
|--|---------------------------------|---|
| Model | : SST-UL | |
| Ventilation system | : RCRv75/2m | |
| Flow | : 10 to 25 m³/h | |
| Total pressure | : Max. 250 Pa | |
| Sound level | : 38 dB(A) at a distance of 3 m | |
| Monitored parameters | : - Flow | <div style="border: 1px solid black; padding: 5px; width: fit-content;"> <p style="text-align: center; margin: 0;">Test</p> <p style="text-align: center; margin: 0;">Volumenstrom Flow</p> <div style="display: flex; justify-content: space-around; margin: 0;"> ○○○ </div> <p style="text-align: center; margin: 0;">NC COMM NO</p> <p style="text-align: center; margin: 0;">Sättigung Saturation</p> <div style="display: flex; justify-content: space-around; margin: 0;"> ○○○ </div> <p style="text-align: center; margin: 0;">NC COMM NO</p> <p style="text-align: center; margin: 0;">230V/50Hz</p> <p style="text-align: center; margin: 0;">In=0,2A</p> <p style="text-align: center; margin: 0;">Fuse0,6A</p> </div> |
| - Saturation of activated carbon filter | | |
| LEDs | : ROT = Error, GREEN = Service | |
| | : 25 / 50 / 75% | |
| | : Ambient concentration | |
| Signal outputs | : - Potential-free for flow | |
| - Potential-free for saturation of activated carbon filter | | |
| Contact power | : 230 V / 8 A | |
| Power supply | : 230 V / 50 Hz | |
| Current consumption | : Max. 0.5 A | |
| IP protection | : CE II _/3G c e IIB T4 | |
| Frequency | : 50 Hz | |
| Motor protection | : PTO 90°C | |
| Weight | : 20 kg | |


Characteristic

Dimensions

Accessories shipped with the SST-UL

Installing the SST-UL

- Insert the adapter
- Place the device on top of the cabinet with the adapter over the exhaust vent, then place the height adjuster underneath the device and align
- Fasten the height adjuster to the cabinet using the screws provided

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